

Our Products

DemeTECH manufactures a full range of high quality surgical and dental sutures. Every batch of sutures passes rigorous quality assurance testing to guarantee our customers the highest grade product.

Absorbable

- Polyglycolic Acid
- Polyglactin 910
- Rapid Absorbable

Non Absorbable

- Nylon Monofilament
- Nylon Multifilament
- Silk
- Polyester braided
- Polyester braided without covering
- Polypropylene




Contact Us

DemeTech Corporation

3530 NW 115 Ave
Miami, FL 33178

Toll Free (888) 324-2447
Phone (305) 597-5277
Fax (305) 437-7607

www.demetech.us
info@demetech.us



Surgical Sutures

About DemeTech

Surgical Sutures

Company Info

DemeTECH Corporation is a world-renowned leader in surgical sutures and blades. DemeTECH strives to enhance doctor-patient relationships through hard work, dedication, and a commitment to the pursuance of cutting-edge technology and innovation. Our corporate headquarters are located in Miami, Florida, with office facilities stemming throughout North America, South America, Asia, the Middle East, and Europe.



DemeTECH is an exceptional corporation, in that it uses next-generation technology to meet the growing demands of today's society. Our quality has been proven time and time again: high demand and use worldwide, registration in over 30 countries, and certification of approval of the three major medical authorities-FDA, CE, and ISO. Through its incredible commitment to efficiency, DemeTECH is able to provide sutures of paramount quality at a lower and more reasonable cost. Thus, with a prosperous future ahead, DemeTECH sees no boundaries to its success and wishes to include you in this flourishing voyage.

Polyglycolic Acid



DemeTECH's Polyglycolic Acid suture is a synthetic, absorbable, sterile, surgical suture composed of 100% glycolide and coated with a blend of polycaprolactone and calcium stearate. DemeTECH's PGA sutures are available dyed and colorless.

Nylon Monofilament




DemeTECH's Nylon suture is a non absorbable, sterile surgical suture composed of the long chain aliphatic polymers Nylon 6 and Nylon 6.6. DemeTECH's Nylon sutures are available in either black or blue, dyed with FDA listed dye(s).

Silk - Braided



DemeTECH's Silk suture is a non-absorbable, sterile, surgical suture composed of an organic protein called fibroin. DemeTECH's silk sutures are processed to remove the natural waxes and gums from its origin. DemeTECH's silk sutures are indicated for use in general soft tissue approximation and/or ligation.

In vivo Tensile Strength (USP 2-0)



* Dexon is a registered trademark of United States Surgical, a division of Tyco Healthcare Group LP. All rights reserved.

† Vicryl is a registered trademark of Ethicon, Inc.

DemeTECH Corporation
3530 NW 115 Ave
Miami, FL 33178

Toll Free (888) 324-2447
Phone (305) 597-5277
Fax (305) 437-7607